

▶ ENDO jars

- 1| Made of polypropylene.
- 2| With spout to avoid spillage and handle to ensure an easy grip.
- 3| With permanent embossed graduation.

- 4| Highly transparent, especially when containing liquid.
- 5| Autoclavable

▶ Graduated jug tall form, PP



Code	Cap.	Grad.	d.	h.
CCB001	250 mL	10/1	82 mm	100 mm
CCB002	500 mL	10/1	95 mm	125 mm
CCB003	1000 mL	100/1	110 mm	180 mm
CCB004	2000 mL	20/1	155 mm	220 mm



▶ Graduated jug short form, PP



Code	Cap.	Grad.	d.	h.
CCB005	250 mL	10/1	85 mm	85 mm
CCB006	500 mL	10/1	90 mm	118 mm
CCB007	1000 mL	50/1	119 mm	130 mm
CCB008	2000 mL	20/1	143mm	167 mm
CCB009	3000 mL	100/1	167 mm	180 mm
CCB010	5000 mL	250/1	190 mm	270 mm

