



High form PP measuring cylinder

Code. DBJ001, DBJ002, DBJ003, DBJ004,
DBJ005, DBJ006, DBJ007, DBJ008

- Made of polypropylene
- With hexagonal base that allows a good stability of the cylinders and all of them present spout for an easy pouring
- Permanently moulded graduations



Technical specifications

Code	Capacity	Inner \varnothing^*	Height*	Divisions
DBJ001	10 mL	12,5 mm	147 mm	0,2 mL
DBJ002	25 mL	18 mm	171 mm	0,5 mL
DBJ003	50 mL	23,5 mm	250 mm	1 mL
DBJ004	100 mL	28 mm	250 mm	1 mL
DBJ005	250 mL	37 mm	328 mm	2 mL
DBJ006	500 mL	52,5 mm	365 mm	5 mL
DBJ007	1000 mL	62 mm	440 mm	10 mL
DBJ008	2000 mL	78,5 mm	536 mm	20 mL

*Approx. dimensions