



▶ PORTABLE METERS - 850 Series

The 850 Series instruments are designed at laboratory level for measurement or analysis of parameters in the work area or field. They combine intelligent functions with a robust structure. Their IP57 protection rating, waterproofing, silicone connection protectors and rubber coating make them suitable for use in harsh environments, such as industrial environments or outdoor measurements with harsh environmental conditions.

They are equipped with a large LCD display for measured values, calibration indicators, stable reading icons and self-diagnostics.

- Parameters: pH, ORP, Conductivity, TDS, Salinity, Temperature
- Supplied with pH and/or conductivity combination electrodes with builtin temperature probe
- Automatic temperature compensation
- Large LCD screen with automatic lock and shutdown
- Battery indicator icons.
- Water and dust resistant IP57
- Self-diagnosis and calibration warning
- Supplied with carrying case, user manual, AA batteries (x3), calibration bottles, 50 mL calibration solutions: pH 7, pH 4 and/or 1413 μ S, 12.88 mS (50 mL)

REF. - KDD007: pH meter

REF. - KZD007: pH meter and conductivity meter

REF. - KDG004: Conductivity meter





▶ PORTABLE METERS - 850 Series - Technical features

Code	KDD007	KZD007	KGD004
Model	pH850	PC850	EC850
Parameter	pH/mV/Temp	pH/mV/Temp/Cond/TDS	Cond/TDS/Temp
pH			
Range	0-14	-	-
Resolution	0.1/0.01 pH	-	-
Accuracy	±0.01 pH ± 1 digit	-	-
Temp. compensation	0-100°C	-	-
Calibration	1 to 3 points (USA/NIST)	-	-
mV			
Range	±1000 mV	-	-
Resolution	1 mV	-	-
Accuracy	± 0.2% FS ± 1 digit	-	-
Conductivity			
Range		0 – 19.99 µS/cm 20.0 – 199.9 µS/cm 200-1999 µS/cm 2 – 19.99 mS/cm 20 - 200.0 mS/cm	
Resolution		0.01/0.1/1 µS, 0.01/0.1 mS	
Accuracy		± 1.0% FS ± 1 digit	
Electrode constant		0.1/1.0/10.0 cm-1	
Temp. compensation		0 – 50°C	
Calibration		1 to 3 points	
TDS			
Range		0.1 mg/mL – 100 g/L	
Coeficiente TDS		0.40-1.00	
Temperature			
Range	0 – 100°C		
Resolution	0.1°C		
Accuracy	±0.5°C ±1 digit		

▶ PORTABLE METERS - 850 Series - Electrodes and conductivity cells



KDG021

Electrode for pH with temperature probe, 201T-F



KZD026

Conductivity electrode with T^a probe (k=1), 2301T-S

Code portable meter	KDD007	KZD007	KDG004
Electrode or conductivity cell	KDG021	KDG021 KZD026	KZD026