


Rotary vane vacuum pump

Code. PXB012

- The NAHITA PXB012 vane vacuum pump is a robust, solid pump. Made of metal with rubber feet and base to avoid slipping due to vibrations. It includes a handle for easy carrying. With general ON/OFF button.
- Oil tank capacity: 330 mL
- Uses vacuum pump oil type HFV-32.



INCLUDES

- Filter for air outlet
- Vacuum pump oil
- Power cable

Technical specifications

Code	PXB012	
Voltage	110V / 60Hz	220V / 50 Hz
Free air displacement	5 CFM	4,5 CFM
Full vacuum	3×10^{-1} Pa 0,05 Mbar	
Rotation speed	3500	2800
Power	1/2HP	
Oil capacity	330 mL	
Dimensions	290x115x220 mm	
Weight	9 Kg	