




Digital, RSLAB-13 PRO

Code LLB003, LLB005

- Anti-spill samples: prevent samples overflow.
- Overload protection: motor would automatically stops in case of continuous overload.
- Real-time display actual speed and setting speed value.
- Extended application: PC control via RS232 interface
- Universal stand (CODE:LLF001) and stirrer (CODE: LLF002)



Technical specifications

Code	LLB003	LLB005
Screen	LCD	LCD
Speed range	50-2200 rpm	50-2200 rpm
Max. volume	20 L	40 L
Max. viscosity	10000 mPas	50000 mPas
Rod dimensions	13 mm	13 mm
Approx. dimensions	83x220x186 mm	83x220x186 mm
Approx weight	2,6Kg	2,8Kg



Digital, RSLAB-13 PRO

ACCESSORIES

Accessories included with the equipment



Soporte universal
Code LLF001*
Complete universal stand
* included with the equipment.



Crossed stirrer
Code LLF002*
5x40 cm, stainless steel
* included with the equipment.

Compatible accessories for the equipment



Straight stirrer
Code LLF003
6x40 cm, stainless steel



Blade stirrer
Code LLF004
7x40 cm, stainless steel



Centrifugal stirrer
Code LLF005
9x40 cm, stainless steel



Crossed stirrer
Code LLF006
6.5x35 cm, PTFE



Straight stirrer
Code LLF007
7x35 cm, PTFE



Blade stirrer
Code LLF008
7x35 cm, PTFE