

Embudos de decantación | Separating funnels

Marca|Trademark: Nahita

Los embudos de decantación se utilizan para la separación de líquidos inmiscibles de distinta densidad, generalmente agua y diversos disolventes orgánicos (éter, cloroformo o alcoholes). Fabricados en vidrio borosilicato resistente al choque térmico, los embudos de decantación de Nahita presentan forma de pera y tapón de vidrio con esmerilado. La llave de paso fabricada en teflón es resistente a la acción de ácidos y bases.

Separating funnels are used to separate immiscible liquids of different densities, such as water and several organic solvents (ether, chloroform or alcohols). Made of thermal shock resistant borosilicate glass, Nahita separating funnels are pear-shaped and present glass ground stopper. The stopcock is made of teflon resistant to acids and bases.



Referencia Code	Capacidad Capacity	Ø Embudo Funnel Ø	Ø Rama Stem Ø	Altura Height	Esmerilado Ground joint
12394210	100 mL	58 mm	8 mm	210 mm	21/23
12394220	250 mL	80 mm	10 mm	300 mm	29/32
12394250	500 mL	100 mm	15 mm	350 mm	29/32