

▶ “ECO” ISO bottles

1| ISO glass bottles “ECO” are made of transparent soda glass and present white graduation and marking area and DIN GL45 thread.
2| The pouring ring, to avoid drops formation, and the cap are made of blue PP (resistant up to 140° C).

3| They also present a studded base for a good stability and abrasion resistance to increase their durability



CODE: CFB001



CODE: CFB003



CODE: CFB001 from CFB004

| Code | Capacity | Ø Base | Height |
|--------|----------|--------|--------|
| CFB001 | 100 ml | 56 mm | 105 mm |
| CFB002 | 250 ml | 70 mm | 143 mm |
| CFB003 | 500 ml | 86 mm | 182 mm |
| CFB004 | 1000 ml | 101 mm | 203 mm |