

## Placas de porcelana con cavidades | Porcelain spot plates

Marca|Trademark: Nahita

Las placas con cavidades de Nahita están fabricadas en porcelana y presentan gran resistencia a la acción de ácidos y bases y al choque térmico. Gracias a su elevada resistencia química son ideales para llevar a cabo pequeñas reacciones químicas mientras que debido a su color blanco son también idóneas para aquellas reacciones en las que se produzcan cambios de color, e incluso para la tinción de pequeños especímenes. La superficie y los laterales de las placas están esmaltados para facilitar su limpieza y los pocillos presentan aproximadamente 20 mm diámetro y 5 mm profundidad.

Nahita spot plates are made of porcelain and are highly resistant to acids and bases and thermal shock. Thanks to their great chemical resistance are ideal to carry out small chemical reactions while due to their white colour are also perfect for colour reactions or even for staining small specimens. Plate surface and sides are glazed to ease cleaning and each well is approximately 20 mm diameter and 5 mm depth.



Referencia Code	Nº cavidades No. wells	Largo Length	Ancho Width
20255060	6	85 mm	54 mm
20255120	12	115 mm	90 mm