

## ▶ Shower and eye washer



- 1| Shower to be use in case of emergency that is combined with the eyewash in a pedestal.
- 2| It is made in stainless resistant to corrosion.
- 3| The shower head as a big diameter, it gives a big volume of water and a circular water curtain.
- 4| The eyewash has two sprays with a big volume of water with low pressure that can bath the eyes and the face of the person that has been injured.
- 5| The shower operation is by a rigid handle and the eye wash has a lateral manual opening that works pressing the lever.

Code	Material
RDB004	Stainless steel 304