



► Stirrer 3 dimensions RSLAB-14

1 | Provides a 3D gyratory motion, ideal for low foaming agitation, DNA extraction, electrophoresis gel staining and decolorization, samples washes, molecular hybridization and immunosedimentation, etc.
2 | 2 models: the first one is analogical with simple plate, the other one is digital with plate and 8 elastic bands.

3 | Compact space-saving design fits easily in the incubator and cryogenic box.
4 | Aluminum-alloyed platform provides strong bearing capability, with non-slip mat holds vessels prevent slipping.
5 | Protection class according to DIN EN60529: IP21.

► Analogic, RSLAB-14 A



1 | Platform with anti-slip coating included. Includes 3 small and 3 large elastic bands

Code	LJL001
Type of movement	Orbital 3D
Shaking angle	7°
Speed range	0 – 80 rpm
Operation type	Continuous
Max. weight of shaking	3 Kg
Platform material	Coated aluminium alloy
Aprox. dimensions	330 x 270 x 130 mm
Aprox. weight.	3,2 Kg

► Digital, RSLAB-14 D




1 | Supplied with platform with 8 elastic bands. Platform with anti-slip coating included.

Code	LJL002
Type of movement	Orbital 3D
Shaking angle	9°
Speed range	10 – 70 rpm
Operation type	Continuous or scheduled
Timer	1 - 1199 min
Max. weight of shaking	5 Kg
Platform material	Coated aluminium alloy
Aprox. dimensions	430 x 360 x 106 mm
Aprox. weight.	9,0 Kg



► Accessories (not included)

			MODELS	
	Code	Description	14A	14D
	LJS006	Universal bar for plate with 3 bars	●	
	LLJS007	Universal bar for plate with 4 bars		●
	LJS008	Plate with 3 bars	●	
	LJS009	Erlenmeyer fixing clip attachment	●	
	LJS010	Plate with 4 bars		●
	LJS013	Anti-slip platform. Dimensions: 24x21,5 cm		
	LJS011	Erlenmeyer fixing clip attachment		●
	LJS012	Platform with 8 rubbers. Dimensions: 30x28,5 cm		○
	LJS001	Fixing clip for Erlenmeyer flask 28 units of 25 mL 9 units of 25 mL	●	●
	LJS002	Fixing clip for Erlenmeyer flask 16 units of 50 mL 9 units of 50 mL	●	●
	LJS003	Fixing clip for Erlenmeyer flask 16 units of 100 mL 9 units of 100 mL	●	●
	LJS004	Fixing clip for Erlenmeyer flask 9* units of 200/250 mL 4 units of 200/250 mL	●	●
	LJS005	Fixing clip for Erlenmeyer flask 9* units of 500 mL 2 units of 500 mL	●	●
* with 3 fixing points instead of 4				

○ Accessory **included** with the indicated equipment

● Accessory **compatible** with the indicated equipment