



Round bottom flask



▶ 1 neck



- 1 | Round bottom spherical flask with frosted mouth.
- 2 | Complies with USP and ISO/DIN 4797 standards except for articles (*).
- 3 | For other volumes or frosted, please consult.

Code	Cap.	Frosted	h.	d.
PBL001	10 mL	14/23	35 mm	70 mm
PBL002	25 mL	14/23	85 mm	41 mm
PBL003*	25 mL	19/26	85 mm	41 mm
PBL004	50 mL	14/23	90 mm	51 mm
PBL005	50 mL	19/26	90 mm	51 mm
PBL006	50 mL	24/29	115 mm	51 mm
PBL007	50 mL	29/32	115 mm	51 mm
PBL008	100 mL	14/23	105 mm	64 mm
PBL009	100 mL	19/26	105 mm	64 mm
PBL010	100 mL	24/29	105 mm	64 mm
PBL011	100 mL	29/32	105 mm	64 mm
PBL036	150 mL	24/29	137 mm	75 mm
PBL037	150 mL	29/32	137 mm	75 mm
PBL012	250 mL	14/23	140 mm	85 mm
PBL013	250 mL	19/26	140 mm	85 mm
PBL014	250 mL	24/29	140 mm	85 mm
PBL015	250 mL	29/32	140 mm	85 mm
PBL038	250 mL	34/35	140 mm	85 mm
PBL016*	500 mL	19/26	163 mm	105 mm
PBL017	500 mL	24/29	163 mm	105 mm
PBL018	500 mL	29/32	163 mm	105 mm
PBL019	1000 mL	24/29	200 mm	131 mm
PBL020	1000 mL	29/32	200 mm	131 mm
PBL041	5000 mL	45/40	305 mm	223 mm
PBL021	2000 mL	24/29	240 mm	166 mm
PBL022	2000 mL	29/32	240 mm	166 mm
PBL032*	2000 mL	34/35	240 mm	165 mm
PBL043	3000 mL	29/32	260 mm	185 mm
PBL040	1000 mL	19/26	200 mm	131 mm
PBL042	10000 mL	45/40	380 mm	279 mm