



- 1 | User-friendly and intuitive moisture analyzer with 5-inch high-definition touch screen.

2 | Silent double action opening. It has an RS232 port to connect a printer. Possibility to print up to 9 groups of data records stored in memory.
- 3 | Equipped with a level gauge and two adjustable leveling feet to compensate for slight unevenness on the operating table.

Halogen Moisture Analyzer



|                          |            |
|--------------------------|------------|
| Code                     | KCD001     |
| Capacity                 | 100g/1mg   |
| Readability              | 50g/5mg    |
| Repeatability            | 1mg        |
| Operation Temperature    | 5°C-40°C   |
| Pan size                 | 100mm      |
| Heating Element          | Halogen    |
| Temp.Range               | 50°C-180°C |
| Temp.Accuracy            | 1°C        |
| Moisture Range           | 0%-100%    |
| Moisture Readability     | 0.01%      |
| Dry weighing range       | 100%-0%    |
| Dry weighing Readability | 0.01%      |