



Tube roller, RSLAB-10



Features

- It allows shaking several types of tubes with different liquid, viscous samples or solids suspension that require minimum airmg.
- 1 | 6 rollers.
 - 2 | Rocking and rolling action for complete mixing.
 - 3 | Analogical regulation speed.
 - 4 | Continuous operation mode.
 - 5 | Protection class according to DIN EN 60529: IP 21

Technical features

Code	LJD002
Model	RSLAB-10
Shaking movement	Rotary and rolling
Speed range	0-70 rpm
Max. shaking weight	5 Kg
Run type	Time/continuous
Dimensions (LxWxH)	530x260x120 mm
Weight	5.1 Kg

MULTIDIMENSIONAL SHAKERS
[L] STIRRING

Note: manufacturer reserves the right to modify features and photographs of the goods.