

Pipeta para vinos | [Wine pipette](#)

La pipeta para vinos es un instrumento de vidrio utilizado en las bodegas como medio de obtención de muestras de vino de una bodega para su posterior cata o análisis. Consiste en un largo tubo de vidrio con dos asas en un extremo para poder sujetar la pipeta cómodamente mientras es introducida en la bodega. En el extremo opuesto a la punta de la pipeta, el tubo está cerrado a excepción de un pequeño orificio sobre el que se coloca el pulgar para evitar que la pipeta se vacíe al sacarla de la bodega.

Marca|**Trademark:** Nahita

Wine pipette is a glass instrument used in wineries to obtain samples of wine directly from the barrel for a later tasting or analysis. It consists of a long glass tube with a pair of handles at one end to hold the pipette while it is introduced into the barrel. At the opposite end of the pipette tip, the tube is closed except for a little hole over which the thumb is placed to avoid the pipette to empty when it is taken out the barrel.



Referencia	Capacidad	Diámetro	Longitud
Code	Capacity	Diameter	Length
11590000	350 mL	30 mm	640 mm