



▶ BDB - With stopcock

- 1 | Endo Glassware offers a complete series of glass desiccators, with or without connection to vacuum pump, to keep samples and reagents free of moisture.
- 2 | Desiccators present flat flange and the lower portion modelled to hold

the porcelain disk over the desiccant agent, and are supplied together with the lid and the porcelain disk.

- 3 | These components are also available individually.

▶ Desiccators with stopcock lid



- 1 | Lid with ground joint and vacuum stopcock.
- 2 | These desiccators can support a vacuum pressure of 6 KPa.

Code	Capacity	Ø Outer	Ø Inner	Height
BDB002	2 L	210 mm	150 mm	220 mm
BDB003	4 L	270 mm	200 mm	300 mm
BDB004	7 L	305 mm	250 mm	360 mm
BDB005	14 L	365 mm	300 mm	420 mm

▶ BDD - Without stopcock

▶ Desiccators with knob lid




- 1 | Without connection to vacuum pump.
- 2 | The lid presents a big aniel for an easy manipulation.

Code	Capacity	Ø Outer	Ø Inner	Height
BDD001	1 L	136 mm	100 mm	180 mm
BDD002	2 L	210 mm	150 mm	220 mm
BDD003	4 L	270 mm	200 mm	300 mm
BDD004	7 L	305 mm	250 mm	360 mm
BDD005	14 L	365 mm	300 mm	420 mm

▶ BDG - Disks

▶ Porcelain disks for desiccators



- 1 | Nahita brand, glazed on one side.

Code	Diameter	Ø Hole
BDG001	90 mm	10 mm
BDG002	140 mm	10 mm
BDG003	190 mm	10 mm
BDG004	220 mm	10 mm
BDG005	275 mm	10 mm



BDL - Accesorios

▶ Knob lids for desiccators

Code	For desiccators of	Ø Outer	Ø Inner
BDL001	100 mm	136 mm	100 mm
BDL002	150 mm	210 mm	150 mm
BDL003	200 mm	270 mm	210 mm
BDL004	250 mm	305 mm	250 mm
BDL005	300 mm	365 mm	305 mm



▶ Ground stopcock for desiccators

1 | With cone.

Code	Ground
BDL011	24/29
BDL012	19/28



▶ Ground socket lids for desiccators

1 | With socket.

Code	For desiccators of	Ground
BDL013	100 mm	19/28
BDL014	150 mm	19/28
BDL015	200 mm	24/29
BDL016	250 mm	24/29
BDL017	300 mm	24/29

