



Ultrasonic baths

Code ZCC016, ZCC017, ZCC018, ZCC019, ZCC020, ZCC021, ZCC022

- Digital display with temperature and time selection.
- Includes lid and basket with handles for a better grip.
- AISI304 stainless steel water resistant tank.
- Adjustable timer from 1 to 99 min. Selectable from minute to minute
- Adjustable temperature from room to 80 °C
- 220V power supply



Technical specifications

Code	ZCC016	ZCC017	ZCC019	ZCC020	ZCC021	ZCC018	ZCC022
Frequency (KHz)	40	40	40	40	40	40	40
Capacity (L)	2	3	6	10	15	22	30
Tank size (mm)	150x140x100	240x140x100	300x155x150	300x240x150	330x300x150	500x300x150	500x300x200
Unit size (mm)	190x170x195	270x160x220	330x180x280	330x270x280	360x330x280	530x330x280	530x330x330
Heating power (W)	100	100	300	300	400	500	500
Ultrasonic power (W)	60	120	180	240	360	480	600