



Power supply for electrophoresis cell

Code ZFD014, ZFD016

- Easy to use the LCD display and control panel
- Constant voltage, current or power.
- Programed value and actual value are indicated at the same time
- Programmable time function
- Memory for 10 electrophoresis programs
- 4 parallel outputs to connect up to 4 cuvettes at the same time.
- With automatic detection of lack of load, overload, short circuit, rapid change of resistance and system overheating



Technical specifications

Code	ZFD014	ZFD016
Output voltaje	10-300 V	10 - 600 V
Output current	1-400 mA	1 - 500 mA
Output power	1-120 W	1 - 300 W
Resolution	1 V, 1 mA, 1 W	
Timer	1 min to 99 h 59 min	
Backlit	LCD display	
Outputs	4 parallel outputs	
Dimensions	230x237x118 mm	
Power supply	220 V; 50 Hz	