



Uv transiluminator - Main features

- 1 | Uv illumination source wavelength 312 nm
- 2 | Equipped with 4 quartz lamps of 8 W of long life and that allow a uniform illumination.
- 3 | On/off switch on the front of the unit.
- 4 | UV glass filter of 150 x 200mm
- 5 | Screen highly transparent and protective against UV emission. It can be positioned at any angle for safe and comfortable use of equipment.



Code	ZFD015
Wavelength	312 nm
Filter size, UV	150x200 mm
Emission power	4 x 8 W
Dimensions (LxWxH)	325x265x150 mm
Power supply	220 V; 50 Hz